**felisa**<sup>®</sup>

FABRICANTES FELIGNEO, S.A. DE C.V.  
Alfonso Garzón Santibañez No. 7 , Col. Indígena San Juan de Ocotán  
C.P. 45019, Zapopan, Jalisco, México.  
Tel. 01 33) 31 10 60 02, 31 10 60 77, Fax. 31 10 61 03  
[Http://www.felisa.com.mx](http://www.felisa.com.mx)

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
## DIAGRAMA ELECTRICO

# OWNER'S MANUAL FOR LOW TEMPERATURE INCUBATORS

The care you take in reading and following this instruction will probably determinate the satisfactory service you will receive from your Low Temperature Incubator.

## UNPACKING

Carefully remove the Incubator from the shipping case. Preserve all paper work for future reference. If damage has occurred from shipment a claim must be filled with the carrier immediately; preserve the shipping container for inspection by the carrier. Contact your dealer or Felisa.

## INTRODUCTION

FELISA Low Temperature Incubators are equipped with forced air circulation and with a digital microprocessor-based control with high quality brilliancy display and with a resolution of 0.1 ° C. You can operate the Incubator with high and low temperature limit alarms. . All models operates with a electronic sensor. A model with programable light inside the chamber is available.

## INSTALLATION

Install your Incubator using a connection of the adequate voltage. See the label specification before connecting. The Incubator can be installed on any surface or table sufficiently firm and strong, leaving a minimum space of 15 cm between the Incubator and any other vertical surface. Keep the area around free of any material to allow ventilation in the area of the bottom. For correct operation it is necessary that the user becomes familiarized with all the available controls and specifications shown in each model.

## LOAD IN THE CHAMBER

### OPERATION TEMPERATURE

Your Incubator is designed to operate satisfactorily to temperatures that do not exceed 50 ° C, If this temperature is not exceeded the Incubator will provide a long trustable service. The minimum operation temperature in the Incubator can't be less than 2 °C.

## OPERATION

### a) Auto-revision and turn on of the system

The Felisa Low Temperature Incubators has an error auto-detection that let you know from the most common failures of the system, therefore the user can be sure that the Incubator is working in optimal conditions. Every time you turn on the Incubator, the system auto test most of its components. This test is a revision of the crucial features of the system.

The display present for 2 seconds the number "8888" at the same time, the audible alarm is activated. After this, the display turns off and the system starts the internal test. If all the features are all right then the Incubator will present in the display the chamber temperature. If any of the features fail, the

Failure Code:	
Sensor Failure	Err (0 y 4)
Heating Element Failure	Err (2)
Memory Failure	Err (1 y 7)
Overheating	Err (3)

### **b) Setting the working temperature**

The Incubator has a simple menu that allows you to configure the working mode of the system. This menu is accessible when the key switch (27) is in the open position, otherwise you won't be able to change the configuration that is established at that moment.

Press the menu button, the display will present the message "SP" (Set Point) for two seconds and then configure the temperature. With the increase or decrease buttons you can modify this parameter to the desired operation temperature, then press the menu button to continue with the configuration and to record the selected value.

### **c) Setting the overshoot alarm**

At this point the display will present the message "SPHA" and then will present the temperature of the overshoot alarm. This temperature will always be 0.5 ° C higher than the set point temperature.

To modify the overshoot alarm use the increase or decrease buttons to set the desired temperature, and then press the menu button to record and continue with the configuration.

### **d) Setting the over cool alarm**

When you enter at this point the display will present the message "SPAL" and then the temperature of the overcool alarm, this temperature always be 0.5 ° C lower than the set point temperature. Use the increase or decrease buttons to set the desired temperature, and then press the menu button to continue with the configuration.

### **e) Setting the Audible Alarm**

The Incubator has an audible alarm that can be activated if you wish, this alarm will sound when a failure is presented in the system, as well as the display present the error message.

When you enter this section in the menu, the display will present the message "BEEP" for two seconds and then the present status of the alarm ("On" or "Off"). With the increase or decrease buttons you can select the configuration. When you finish your selection press the menu button, the display will turn off, will be sound a little beep indicating that the configuration is finished.

## **OPERATION SEQUENCE**

The system needs some time to reach the operation temperature and then stabilize. Approximately twenty five minutes after the system reaches the set point temperature, the Ok led will turn ON, this indicates that all the alarms are activated and the timer will start working in case you'd configured it previously. When the time is finished the audible alarm will sound, and the system will continue controlling the heating. To turn off the alarm, press the decrease button.

If the Ok led is ON and the temperature surpasses the overshoot or over cool temperature, the system will present the error message, the alarm will turn on and the Incubator will stop the temperature controlling.

## CALIBRATION

Your Incubator was calibrated in the plant with certificated references. However, transportation, ambient temperature or specific operation conditions might affect these calibration. If you have a reliable thermometer you can calibrate your Incubator.

Program the Incubator at your operation temperature, put the thermometer inside the chamber, wait until the temperature is stable and the Ok led is activated. At this point, check the thermometer temperature against the display, if there is a difference you have to calibrate your Incubator. Turn off the Incubator with the on-off switch (29), with the menu button pressed turn on the Incubator, the display will show "CAL" and then a temperature, adjust this number to the thermometer reading using the increase or decrease buttons and then press the menu button. Now your Incubator is calibrated with your thermometer reference.

## 24 HOUR DIGITAL TIME SWITCH

### KEY FUNCTIONS:

DAY. Press to select the right day of the week in mode CLK (clock)

HOUR. Press to put the hour in mode CLK and mode SCH (schedule)

MIN. Press to put the minutes in mode CLK and mode SCH.

CONTROL. Press to change the actual load state (works in mode AUTO and MANUAL).

MODE. Press to change to next mode. The sequence is: CLK, SCH, AUTO, MAN. In modes CLK and SCH the equipment returns automatically to mode AUTO if no enter is given during 5 minutes.

DELETE. Press to delete the information presented when you are in mode AUTO.

AUTO. Press to change from any mode to mode AUTO.

ENTER. Press to storage information presented in memory. The information won't be stored until press the ENTER key.

NOTE: During the programming, each time you press the key you move forward one number. To move faster keep the key pressed.

### SETTING THE CLOCK

Press mode key until CLK appears. Press DAY key to select the present day, press HOUR key to select the present hour. Check the existence of AM or PM. Press MIN key to select the present minute. Press ENTER key and the clock is ready in time and present date.

### SETTING THE PROGRAM

Press MODE key until CLK appears. There are 7 on and 7 off operations, one followed by the other, ( the first is On, the second Off, the third On, etc.) The display will indicate it is ready to receive SCH set point #1 which is an On event. Press HOUR and MIN keys for the first On setting (check for AM or PM).

Press ENTER key to store the information and display will indicate it is ready to receive SCH set point # 2 which is an Off event.

Press HOUR and MIN keys for the first Off setting. Then press ENTER key to store the information and display will indicate it is ready to receive SCH set point # 3 which is the next On event. Proceed for up to 7 On and 7 Off set points. Then press MODE key and unit will go to the AUTO run mode. In the AUTO

## REVIEW / MODIFY

### CLOCK

Press mode key until display shows CLK. Press Day, Hour and Min keys to change to the correct time. Press Enter.

### SCHEDULE

Press Mode key until display shows SCH. Press Enter key repeatedly to review all the scheduled entries. During the review any set point can be modified or deleted. Use the Hour and Min keys to modify. Use the Delete key to eliminate that event. Press Enter key after each modification.

### OVERRIDE

**TEMPORARY:** in order to temporarily change the current On or Off status of a load, simply press the override key when the unit is in the Auto mode. The altered load status will flash. Override is in effect until the next scheduled event. In order to cancel the override, press override key again.

**LONG TERM:** Press Mode key until display shows MAN (manual). Normal schedule will not be executed and the load status will remain unchanged as indicated. Press the Override key to obtain the correct status or to alter the load status.

## REPAIRS

Like any other manufactured product some parts of the Incubator could be damaged after using for a long time. To replace them always use original factory parts, a list of which is included in this User´s Manual. Always order the parts with its corresponding number. All the parts can be ordered from any Distributor or directly from our Plant.

## IMPORTANT

- Do not lean to much the Incubator, and be carefull with the fan in the back lower side of the equipment.
- Before to carry out any maintenance work disconnect the Oven from its source of energy.
- Do not spill any solution inside the camber. Take care with metallic objects inside the upper part of chamber: could be cause of electric shock.
- Do not change the position of the temperature sensor.
- Do not use your Incubator without a ground connection.
- Do not saturate the chamber with material.
- Voltage variations can damage electronic components.
- If the plug cable is damaged a qualified person is needed to change it to avoid damage.

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All products manufactured by Felisa are guaranteed for one year from date of shipment from Felisa plant. All those products returned within one year will be rebuilt or replaced under the guarantee regardless of reason for failure. Electric motors are guaranteed according to the manufacturer policies. Exceptions may be made by Felisa on particular applications where experience has indicated so conditions are so unusual that premature failure can be expected. Transportaion charges in all cases will be at customers expense. Maximum liability in no case will exceed the value of the Felisa product involved.

Felisa has the right to change or modify the products in order to improve their presentation and/or operation

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